




# HOPETOWN

D E N T A L   L A B




# Product Catalog 2018

# Two ways to find what you need!

## 1. Table of Contents (Pg. 3-4)

Need a Rubber diamond polisher?

Just go to the [Table of Contents](#) and click “Rubber Diamond Polishers” to jump to that page! No need to scroll!




ROTARY INSTRUMENTS FOR ZIRCONIA	
<b>Grinders for Zirconia</b>	
CDZY: Yellow Ceramic Diamond Grinders, Coarse	17
CDZ: Green Ceramic Diamond Grinders, Mid-Coarse	17
CDZB: Black Ceramic Diamond Grinders, Medium	17
Diamond Grinding Discs for Zirconia	20
<b>Polishers for Zirconia</b>	
SDZ: Silicone Diamond Polishers	18
RDZ: Rubber Diamond Polishers	18-19
Bristle Brushes	20
Felt Wheels	20

## 2. Search by Pre-Sinter Processing and Post-Sinter Processing

Not sure which instruments are best for a Procedure, such as shaping occlusal anatomy for zirconia?

Check out the [Pre-Sinter Processing and Post-Sinter Processing sections \(Pg. 5-9\)](#) to find out what tools work best for each process to save you time and get the most use out of your Instruments!



**Anatomy**

**Luna Zirconia**

- High Speed / Friction Grip Diamond Burs, **Medium-Coarse**  
171 types available  
Pg. 24-26

**Orea Zirconia**

- High Speed / Friction Grip Diamond Burs, **Medium-Coarse**  
171 types available  
Pg. 24-26
- Low-Speed / Straight Handpiece Lab Burs, **Medium-Coarse**  
66 types available  
Pg. 21-22

# Table of Contents

## Lab Products

FDA / ISO 13485 / ISO 9000 / CE 0197

### MATERIALS

#### Zirconia

Blocks.....	11
Discs (95mm, 98mm, 100mm) .....	12-14

#### PMMA

Blocks .....	15
Discs .....	15

#### Wax

Discs .....	15
-------------	----

Sintering Beads.....	16
Polishing Compound and Diamond Paste .....	20

### ROTARY INSTRUMENTS FOR ZIRCONIA

#### Grinders for Zirconia

CDZY: Yellow Ceramic Diamond Grinders, Coarse.....	17
CDZ: Green Ceramic Diamond Grinders, Mid-Coarse .....	17
CDZB: Black Ceramic Diamond Grinders, Medium .....	17

#### Polishers for Zirconia

SDZ: Silicone Diamond Polishers.....	18
RDZ: Rubber Diamond Polishers.....	18-19
Polishing Brushes, Felt Wheels, Chamois Wheels, Shaping Brushes, Mandrels, Diamond Dresser.....	20

#### Diamonds for Zirconia

Diamond Cutting Discs for Zirconia .....	21
Low Speed Lab Diamond Burs .....	22-23
High Speed / Friction Grip Diamond Burs.....	25-29

## Carbide Burs

Low Speed Lab Carbide Burs .....	24
High Speed / Friction Grip Carbide Burs.....	30

## ROTARY INSTRUMENTS FOR CERAMICS

SPC: Silicone Polishers for Ceramics .....	31-32
--	-------

## ROTARY INSTRUMENTS FOR ALLOYS

SPA: Silicone Polishers for Noble / High Noble Alloys.....	33-34
RPA: Rubber Polishers for Alloys.....	35-37

## ROTARY INSTRUMENTS FOR ACRYLIC AND RESIN

SPR: Silicone Polishers for Acrylic.....	38
--	----

# Clinical Products

FDA / ISO 13485 / ISO 9000 / CE 0197

## CLINICAL ROTARY INSTRUMENTS

Clinical Chairside Kits .....	39-40
-------------------------------	-------

## Clinical Grinders

Ceramic Diamond Grinders for Zirconia (CDZ).....	41
Ceramic Grinders for Ceramic / Resin / Alloy (CCA/ACA).....	41

## Clinical Polishers

Rubber Diamond Polishers for Zirconia & Ceramic .....	42
Rubber Polishers for Alloys .....	43
Silicone Diamond Polishers for Resin.....	45
Silicone Polishers for Ceramic & Resin .....	46
Silicone Polishers for Noble/High Noble Alloy .....	47

## Clinical Diamond and Carbide Burs

High Speed / Friction Grip Diamond Burs.....	25-29
High Speed / Friction Grip Carbide Burs.....	30

## Clinical Equipment

OpaScan: Digital Dental System .....	48
--------------------------------------	----

## Clinical Materials

OpaFit Impression Material .....	49
----------------------------------	----

## Removing Crown from Disc/Block



### For Dry Milling (Discs)



Click to go to  
Diamond  
Burs

High Speed / Friction  
Grip Diamond Burs,  
**Medium-Coarse**  
171 types available  
Pg. 25-27



Click to go to  
Carbide  
Burs

High Speed / Friction Grip  
Carbide Burs  
34 types available  
Pg. 30



### For Wet Milling (Blocks)



Click to go to Diamond Cutting  
Discs

Diamond Cutting Discs  
32 types available  
Pg. 21

## Smooth Sprue Area



Click to go to  
Diamond  
Burs

High Speed / Friction  
Grip Diamond Burs,  
**Medium-Coarse**  
171 types available  
Pg. 25-27



Click to go to  
SPCs

Silicone Polishers  
(SPC)  
69 types available  
Pg. 31-32

## Crown Gross Shaping



Click to go  
to  
Shaping  
Brushes

Shaping Brushes  
(BRW)  
3 types available  
Pg. 20



Click to go to  
SPCs

Silicone Polishers  
(SPC)  
69 types available  
Pg. 31-32

## Occlusal Anatomy



[Click to go to  
Diamond  
Burs](#)

High Speed / Friction  
Grip Diamond Burs,  
**Medium-Coarse**  
171 types available  
Pg. 25-27



[Click to go to  
Carbide  
Burs](#)

High Speed / Friction  
Grip Carbide Burs  
34 types available  
Pg. 30

## Smoothing and Finishing



[Click to go to  
Shaping  
Brushes](#)

Shaping Brushes  
(BRW)  
3 types available  
Pg. 20

## Interproximal Finishing



[Click to go to  
SPCs](#)

Silicone Polishers  
(SPC)  
69 types available  
Pg. 31-32

## Sintering



[Click to go to  
Sintering  
Beads](#)

Zirconia Sintering Beads  
Available in 200gr. or 1kg.  
Pg. 16

## Seat Restoration

### Luna Zirconia



[Click to go to CDZs](#)

**Green** Ceramic  
Diamond Grinders,  
Mid-Coarse (CDZ)  
19 types available  
Pg. 17



[Click to go to Diamond Burs](#)

High Speed / Friction  
Grip Diamond Burs,  
**Medium-Coarse**  
171 types available  
Pg. 25-27



[Click to go to Diamond Burs](#)

Low-Speed / Straight  
Handpiece Lab Burs,  
**Medium-Coarse**  
66 types available  
Pg. 22-23

### Orea Zirconia



[Click to go to CDZs](#)

**Green** Ceramic  
Diamond Grinders,  
Mid-Coarse (CDZ)  
19 types available  
Pg. 17



[Click to go to Diamond Burs](#)

High Speed / Friction  
Grip Diamond Burs,  
**Fine**  
47 types available  
Pg. 28-29



[Click to go to Diamond Burs](#)

High Speed / Friction  
Grip Diamond Burs,  
**Medium-Coarse**  
171 types available  
Pg. 25-27



[Click to go to Diamond Burs](#)

Low-Speed / Straight  
Handpiece Lab Burs,  
**Medium-Coarse**  
66 types available  
Pg. 22-23

## Finish Margins

### Luna Zirconia



[Click to go to CDZs](#)

**Green** Ceramic  
Diamond Grinders,  
Mid-Coarse (CDZ)  
19 types available  
Pg. 17



[Click to go to CDZs](#)

**Yellow** Ceramic  
Diamond Grinders,  
Coarse  
5 types available  
Pg. 17



[Click to go to Diamond Burs](#)

Low-Speed / Straight  
Handpiece Lab Burs,  
**Medium-Coarse**  
66 types available  
Pg. 22-23

### Orea Zirconia



[Click to go to SDZs](#)

Silicone Diamond  
Polishers (SDZ)  
12 types available  
Pg. 18



[Click to go to RDZs](#)

Rubber Diamond  
Polishers (RDZ)  
49 types available  
Pg. 18-19



[Click to go to Diamond Burs](#)

Low-Speed / Straight  
Handpiece Lab Burs,  
**Medium-Coarse**  
66 types available  
Pg. 22-23

## Gross Finishing

### Luna Zirconia



[Click to go to CDZs](#)

**Green** Ceramic  
Diamond Grinders,  
Mid-Coarse (CDZ)  
19 types available  
Pg. 17



[Click to go to CDZs](#)

**Yellow** Ceramic  
Diamond Grinders,  
Coarse  
5 types available  
Pg. 17

### Orea Zirconia



[Click to go to CDZs](#)

**Black** Ceramic  
Diamond Grinders,  
Coarse  
5 types available  
Pg. 17



[Click to go to SDZs](#)

**Silicone** Diamond  
Polishers (SDZ)  
12 types available  
Pg. 18

## Contacts and Occlusions

### Luna Zirconia



[Click to go to CDZs](#)

**Green** Ceramic  
Diamond Grinders,  
Mid-Coarse (CDZ)  
19 types available  
Pg. 17



[Click to go to CDZs](#)

**Yellow** Ceramic  
Diamond Grinders,  
Coarse  
5 types available  
Pg. 17

### Orea Zirconia



[Click to go to CDZs](#)

**Black** Ceramic  
Diamond Grinders,  
Coarse  
5 types available  
Pg. 17



[Click to go to SDZs](#)

**Silicone** Diamond  
Polishers (SDZ)  
12 types available  
Pg. 18



## Anatomy

### Luna Zirconia



[Click to go to Diamond Burs](#)

High Speed / Friction Grip Diamond Burs, **Medium-Coarse**  
171 types available  
Pg. 25-27

### Orea Zirconia



[Click to go to Diamond Burs](#)

High Speed / Friction Grip Diamond Burs, **Medium-Coarse**  
171 types available  
Pg. 25-27



[Click to go to Diamond Burs](#)

Low-Speed / Straight Handpiece Lab Burs, **Medium-Coarse**  
66 types available  
Pg. 22-23

## Bridge Connector Areas

### Luna Zirconia



[Click to go to CDZs](#)

**Green** Ceramic Diamond Grinders, Mid-Coarse (CDZ)  
19 types available  
Pg. 17



[Click to go to CDZs](#)

**Yellow** Ceramic Diamond Grinders, Coarse  
5 types available  
Pg. 17



[Click to go to Sintered Diamond Discs](#)

Sintered Diamond Discs, 7 types  
Pg. 21



[Click to go to Diamond Burs](#)

Low-Speed / Straight Handpiece Lab Burs, **Medium-Coarse**  
66 types available  
Pg. 22-23

### Orea Zirconia



[Click to go to Diamond Burs](#)

Low-Speed / Straight Handpiece Lab Burs, **Medium-Coarse**  
66 types available  
Pg. 22-23



## **Luna™ and Orea™ Zirconia**

Every employee at Opa Dental has worked in laboratories, and all of our products reflect that. We only provide you with the best, high-value products available, along with full technical support. All of our zirconia is designed to our specifications, and has undergone extensive testing before we will put our name on it.

**FDA and CE registered. ISO 13485 Certified.**



### **Luna™ Substructure / Framework Zirconia (1150-1300 MPa.)**

**Recommendations:**


- *Substructure*
- *Substructure full arch bridges*



### **Luna™ Translucent and Pre-Shaded Translucent Zirconia (900-1100 MPa.)**

**Recommendations:**

- *Full contour, Anterior and Posterior*
- *Full arch bridges*



### **Orea Plus™ Multi-Layered Zirconia (900-1100 MPa.)**

**Recommendations:**

- *All Indications*
- *Crowns and bridges*

**ZIRCONIA / BLOCKS**



**Luna™ Substructure / Framework Zirconia Blocks for Sirona Machines (1150-1300 MPa.)**

<b>20x15</b> 12 per Box
<b>20x19</b> 10 per Box
<b>40x15</b> 8 per Box
<b>40x19</b> 6 per Box
<b>55x19</b> 5 per Box
<b>65x25</b> 4 per Box
<b>85x40</b> 2 per Box

**Luna™ Unshaded Translucent and Pre-Shaded Translucent Zirconia Blocks for Sirona Machines (900-1100 MPa.)**

	<b>Unshaded Translucent</b>	<b>Pre-Shaded Translucent</b>
<b>20x15</b> 12 per Box	•	•
<b>20x19</b> 10 per Box	•	•
<b>40x15</b> 8 per Box	•	
<b>40x19</b> 6 per Box	•	•
<b>55x19</b> 5 per Box	•	•
<b>65x25</b> 4 per Box	•	•
<b>85x40</b> 2 per Box	•	•



**Orea Plus™ Multi-Layered Zirconia Blocks for Sirona Machines (900-1100 MPa.)**

	<b>Unshaded (Bleach)</b>	<b>Multilayered Pre-Shaded</b>
<b>20x19</b> 10 per Box	•	•
<b>40x19</b> 6 per Box	•	•

**ZIRCONIA / 95mm Discs**

<b>Luna™ Substructure / Framework 95mm Zirconia Discs (1150-1300 MPa.)</b>	
<b>95x12mm</b>	•
<b>95x14mm</b>	•
<b>95x16mm</b>	•
<b>95x18mm</b>	•

<b>Luna™ Unshaded Translucent and Pre-Shaded Translucent 95mm Zirconia Discs (900-1100 MPa.)</b>		
	<b>Unshaded Translucent</b>	<b>Pre-Shaded Translucent</b>
<b>95x10mm</b>	•	
<b>95x12mm</b>	•	
<b>95x14mm</b>	•	•
<b>95x16mm</b>	•	
<b>95x18mm</b>	•	•
<b>95x20mm</b>	•	•
<b>95x22mm</b>	•	•
<b>95x25mm</b>	•	•
<b>95x30mm</b>	•	•

<b>Luna-Orea™ Cosmetic 95mm Zirconia Discs (650-800 MPa.)</b>		
	<b>Orea Unshaded</b>	<b>Orea Multilayered Pre-Shaded</b>
<b>95x12mm</b>	•	
<b>95x14mm</b>	•	
<b>95x18mm</b>	•	

**ZIRCONIA / 98mm Discs**

**Luna™ Substructure / Framework  
98mm Zirconia Discs  
(1150-1300 MPa.)**

98x10mm	•
98x12mm	•
98x14mm	•
98x16mm	•
98x18mm	•
98x20mm	•
98x25mm	•

**Luna™ Unshaded Translucent and  
Pre-Shaded Translucent  
98mm Zirconia Discs (900-1100 MPa.)**

	Unshaded Translucent	Pre-Shaded Translucent
98x10mm	•	•
98x12mm	•	•
98x14mm	•	•
98x16mm	•	•
98x18mm	•	•
98x20mm	•	•
98x22mm	•	
98x25mm	•	•
98x30mm	•	



**Orea Plus™ Multi-Layered  
98mm Zirconia Discs (900-1100 MPa.)**

	Unshaded (Bleach)	Multi-layered Pre-Shaded
98x12mm	•	•
98x14mm	•	•
98x16mm	•	•
98x18mm	•	•
98x20mm	•	•
98x22mm	•	
98x25mm	•	•

# ZIRCONIA / 100mm Discs

**Luna™ Substructure / Framework  
100mm Zirconia Discs  
(1150-1300 MPa.)**

100x10mm	•
100x12mm	•
100x14mm	•
100x16mm	•
100x18mm	•
100x20mm	•
100x25mm	•

**Luna™ Unshaded Translucent and  
Pre-Shaded Translucent  
100mm Zirconia Discs (900-1100 MPa.)**

	Unshaded Translucent	Pre-Shaded Translucent
100x10mm	•	
100x12mm	•	•
100x14mm	•	•
100x16mm	•	•
100x18mm	•	•
100x20mm	•	
100x22mm	•	•
100x25mm	•	•

**Luna-Orea™ Cosmetic  
100mm Zirconia Discs (650-800 MPa.)**

	Orea Unshaded	Orea Multilayered Pre-Shaded
100x12mm	•	
100x14mm	•	•
100x16mm	•	
100x18mm	•	



# PMMA

## Blocks for Sirona Machines

**40x19**

Clear, A0, A1, A2, A3, A3.5, A4, B1, C2, D2, D3

**55x19**

Clear, A1, A2, A3, A3.5, B1, C2, C3

## 95mm Discs (Zirkonzahn)

**95x16**

Clear, A1, A2, A3, A4, B1, B2, C1, C2, C4, D2, D3

**95x20**

Clear, A1, A2, A3, A4, B1, B2, B3, C1, C2, C4, D2, D3

## 98mm (Roland, Wieland, MCX5, etc.)

**98x14**

Clear, A0, A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4

**98x20**

Clear, A0, A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4

**98x25**

Clear

# WAX



## 95mm Discs (Zirkonzahn)

**95x14**

Blue,  
White

## 98mm (Roland, Wieland, MCX5, etc.)

**98x14**

Blue,

White

**98x20**

White

## 100mm (DENTMILL)

**100x14**

Blue

**ZIRCONIA / Sintering**

**Beads**



**Zirconia Sintering Beads, 2.3mm**






**200 grams**

**1 kilogram**



# ZIRCONIA / Ceramic Diamond Grinders

## Yellow Ceramic Diamond Grinders, Coarse (Recommended for Luna Zirconia)

					
Item	CDZY3C	CDZY4C	CDZY5C	CDZY7C	CDZY11C
Size (MM)	22 x 2	16 x 3.5	22 x 4.5	22 x 2	4 x 10.5
RPM			20,000-35,000		
Price					

## Green Ceramic Diamond Grinders, Mid-Coarse (Recommended for Luna Zirconia)

							
Item	CDZ1MC	CDZ2MC	CDZ3MC	CDZ4MC	CDZ5MC	CDZ6MC	CDZ7MC
Size (MM)	13 x 7	16 x 2	22 x 2	16 x 3.5	22 x 4.5	16 x 1.5	22 x 2
RPM				15,000-30,000			
Price							

							
Item	CDZ8MC	CDZ9MC	CDZ10MC	CDZ11MC	CDZ12MC	CDZ13MC	CDZ14MC
Size (MM)	3 x 7	5 x 12	7 x 12	4 x 10.5	3 x 7	5 x 2.5	13 x 3.5
RPM				20,000-35,000			
Price							

					
Item	CDZ15MC	CDZ16MC	CDZ17MC	CDZ18MC	CDZ20MC
Size (MM)	3.6 x 7	6 x 8	12 x 2	12 x 3	3 x 3
RPM			20,000-35,000		
Price					

## Black Ceramic Diamond Grinders, Medium (Recommended for Orea Zirconia)

					
Item	CDZB1M	CDZB3M	CDZB4M	CDZB5M	CDZB7M
Size (MM)	13 x 7	22 x 2	16 x 3.5	22 x 4.5	22 x 2
RPM			20,000-35,000		
Price					



## Silicone Diamond Polishers for Zirconia (SDZ)

**ZIRCONIA / Silicone  
Diamond Polishers**



Item					
	<b>SDZ1EC</b> Extra Coarse 22 x 3.4	<b>SDZ1C</b> Coarse 22 x 3.4	<b>SDZ1M</b> Medium 22 x 3.4	<b>SDZ1F</b> Fine 22 x 3.4	<b>SDZ1UF</b> Ultra Fine 22 x 3.4
Size (MM)					
RPM			<b>15,000-25,000</b>		
Price					
Item					
	<b>SDZ7EC</b> Extra Coarse 26 x 2.1	<b>SDZ7C</b> Coarse 26 x 2.1	<b>SDZ7M</b> Medium 26 x 2.1	<b>SDZ7F</b> Fine 26 x 2.1	<b>SDZ7UF</b> Ultra Fine 26 x 2.1
Size (MM)					
RPM			<b>15,000-25,000</b>		
Price					

## Rubber Diamond Polishers for Zirconia (RDZ)

**ZIRCONIA / Rubber Diamond  
Polishers**

					
	Extra Coarse	Mid-Coarse	Fine	Extra Fine	Ultra Fine
Item	<b>RDZW1EC</b>	<b>RDZW3MC</b>	<b>RDZW4F</b>	<b>RDZW5EF</b>	<b>RDZ6UF</b>
Size (MM)	26 x 2.1	26 x 2.1	26 x 2.1	26 x 2.1	26 x 2.1
RPM			<b>15,000-25,000</b>		
Price					







### Mid-Coarse Grit

							
Item	<b>RDZ1MC</b>	<b>RDZ2MC</b>	<b>RDZ3MC</b>	<b>RDZ4MC</b>	<b>RDZ5MC</b>	<b>RDZ6MC</b>	<b>RDZ7MC</b>
Size (MM)	3 x 7	5 x 12	3 x 7	4 x 12	4 x 12	9 x 8.5	12 x 1
RPM				<b>15,000-25,000</b>			
Price							
							
Item	<b>RDZ8MC</b>	<b>RDZ9MC</b>	<b>RDZ10MC</b>	<b>RDZ11MC</b>			
Size (MM)	12 x 2	12 x 4	12 x 2	18 x 3			
RPM		<b>15,000-25,000</b>					
Price							










**Rubber Diamond Polishers for Zirconia  
Continued on Next Page**

# ZIRCONIA / RUBBER DIAMOND POLISHERS

## Fine Grit

							
Item	RDZ1F	RDZ2F	RDZ3F	RDZ4F	RDZ5F	RDZ6F	RDZ7F
Size (MM)	3 x 7	5 x 12	3 x 7	4 x 12	4 x 12	9 x 8.5	12 x 1
RPM				15,000-25,000			
Price							
							
Item	RDZ8F	RDZ9F	RDZ10F	RDZ11F			
Size (MM)	12 x 2	12 x 4	12 x 2	18 x 3			
RPM		15,000-25,000					
Price							

## Extra-Fine Grit

							
Item	RDZ1EF	RDZ2EF	RDZ3EF	RDZ4EF	RDZ5EF	RDZ6EF	RDZ7EF
Size (MM)	3 x 7	5 x 12	3 x 7	4 x 12	4 x 12	9 x 8.5	12 x 1
RPM				10,000-20,000			
Price							
							
Item	RDZ8EF	RDZ9EF	RDZ10EF	RDZ11EF			
Size (MM)	12 x 2	12 x 4	12 x 2	18 x 3			
RPM		10,000-20,000					
Price							

## Ultra-Fine Grit

							
Item	RDZ1UF	RDZ2UF	RDZ3UF	RDZ4UF	RDZ5UF	RDZ6UF	RDZ7UF
Size (MM)	3 x 7	5 x 12	3 x 7	4 x 12	4 x 12	9 x 8.5	12 x 1
RPM				10,000-20,000			
Price							
							
Item	RDZ8UF	RDZ9UF	RDZ10UF	RDZ11UF			
Size (MM)	12 x 2	12 x 4	12 x 2	18 x 3			
RPM		10,000-20,000					
Price							

**ZIRCONIA / Bristle Brushes, Felt Wheels, Chamois, Polishing Pastes, Shaping Brushes, Mandrels**

**Polishing Brushes and Wheels**



<b>Item</b>	<b>BRBHMBLK</b> Mounted Stiff Horse Hair Brush	<b>BRBHM</b> Mounted Medium Goat Hair Brush	<b>BRBHMW</b> Mounted Soft, White Goat Hair Brush	<b>FW1</b> Medium Felt Polishing Wheel	<b>CHAMI</b> Mounted Leather Chamois Polishing Wheel
<b>Size (MM)</b>	22	22	22	22 x 6.5	22 x 8
<b>Price</b>					

**Shaping Brushes**



<b>Item</b>	<b>BRW1C</b> Coarse Shaping Brush	<b>BRW1M</b> Medium Shaping Brush	<b>BRW1F</b> Fine Shaping Brush
<b>Size (MM)</b>	25 x 12	25 x 12	25 x 12
<b>Price</b>			

**Polishing Compound and Diamond Paste**



<b>Item</b>	<b>PCHS1</b> High Shine Polishing Compound (Pre-Sinter)	<b>DPP1</b> Diamond High Shine Polishing Paste
<b>Size</b>	55 grams	5 grams
<b>Price</b>		

**Mandrels**



<b>Item</b>	<b>MS1</b> Mandrel, Small Head	<b>MS2</b> Mandrel for Rubber Points	<b>MS3</b> Mandrel, Large Head
<b>Price</b>			




**Diamond Dresser for shaping Grinders and Polishers**






<b>Item</b>	<b>DIDRESS</b> Diamond Dresser for shaping grinders and polishers
<b>Price</b>	

# ZIRCONIA / Diamond Cutting Discs



## Diamond Cutting Discs for Zirconia, 2-Sided

					
	<b>Item</b>	<b>Size (MM)</b>	<b>Item</b>	<b>Size (MM)</b>	<b>Item</b>
	DD1D22.20	22 x 0.20	DD2D22.20	22 x 0.20	DD3D22.20
	DD1D22.25	22 x 0.25	DD2D22.25	22 x 0.25	DD3D22.25
	DD1D22.30	22 x 0.30	DD2D22.30	22 x 0.30	DD3D22.30
	DD1D22.35	22 x 0.35	DD2D22.35	22 x 0.35	DD3D22.35


Price

					
	<b>Item</b>	<b>Size (MM)</b>	<b>Item</b>	<b>Size (MM)</b>	<b>Item</b>
	DD4D22.20	22 x 0.20	DD5D22.20	22 x 0.20	DD6D22.20
	DD4D22.25	22 x 0.25	DD5D22.25	22 x 0.25	DD6D22.25
	DD4D22.30	22 x 0.30	DD5D22.30	22 x 0.30	DD6D22.30
	DD4D22.35	22 x 0.35	DD5D22.35	22 x 0.35	DD6D22.35

Price

			
	<b>Item</b>	<b>Size (MM)</b>	<b>Item</b>
	DD7D22.20	22 x 0.20	DD8D22.20
	DD7D22.25	22 x 0.25	DD8D22.25
	DD7D22.30	22 x 0.30	DD8D22.30
	DD7D22.35	22 x 0.35	DD8D22.35

Price

			
	<b>Item</b>	<b>Size (MM)</b>	<b>Item</b>
	DD35MM	35 x 0.30	DD45MM

Price

## Sintered Diamond Discs for Zirconia, 2-Sided

Coarse Grit

Medium Grit










Fine Grit










					
<b>Item</b>	<b>Item</b>	<b>Size (MM)</b>	<b>Item</b>	<b>Size (MM)</b>	<b>Item</b>
SDDC22.25	SDDM22.25	22 x 0.25	SDDF22.20	22 x 0.20	
SDDC22.30	SDDM22.30	22 x 0.30	SDDF22.25	22 x 0.25	
			SDDF22.30	22 x 0.30	










Price










**ZIRCONIA / Low Speed Diamond Lab Burs**








**Low Speed Diamond Lab Burs  
Mid-Coarse Grit**

									
Item	DBR1 MC	DBR2 MC	DBR3 MC	DBR4 MC	DBR5 MC	DBR6 MC	DBR7 MC	DBR8 MC	DBR9 MC
Size (MM)	1 x 0.7	1.4 x 1.1	1.8 x 1.4	2.3 x 1.9	1.2 x 1.4	1.6 x 1.4	2.1 x 2.2	2.5 x 2.4	2 x 0.8
RPM	20,000-35,000								
Price									

									
Item	DBR10 MC	DBR11 MC	DBR12 MC	DBR13 MC	DBR14 MC	DBR15 MC	DBR16 MC	DBR17 MC	DBR18 MC
Size (MM)	3 x 1.2	4 x 1.4	5 x 2	1.6 x 4	2.1 x 5	2.5 x 6	0.8 x 3	1 x 3.5	1.4 x 3.5
RPM	20,000-35,000								
Price									

									
Item	DBR19 MC	DBR20 MC	DBR21 MC	DBR22 MC	DBR23 MC	DBR24 MC	DBR25 MC	DBR26 MC	DBR27 MC
Size (MM)	2 x 3.5	2.5 x 3.5	3.5 x 3.5	1.2 x 6	1.6 x 6	2.3 x 6	3 x 6	4 x 6	1.2 x 12
RPM	20,000-35,000								
Price									




									
Item	DBR28 MC	MC	DBR30 MC	DBR31 MC	DBR32 MC	DBR33 MC	DBR34 MC	DBR35 MC	DBR36 MC
Size (MM)	1.6 x 12	2.3 x 12	3 x 12	4 x 12	1 x 4	1.4 x 4	1.8 x 4	2.5 x 4	1.2 x 10
RPM	20,000-35,000								
Price									

							
Item	DBR37 MC	DBR38 MC	DBR39 MC	DBR40 MC	DBR41 MC	DBR42 MC	DBR43 MC
Size (MM)	1.6 x 10	2.1 x 10	2.5 x 10	3 x 10	3.5 x 10	4 x 10	1.2 x 8
RPM	20,000-35,000						
Price							








**Low Speed Diamond Lab Burs Continued on Next Page**

ZIRCONIA / Low Speed Diamond Lab Burs




Mid-Coarse Grit

			
Item Size (MM)	DBR44MC 1.8 x 8	DBR45MC 2.5 x 8	DBR46MC 4 x 8


RPM Price 20,000-35,000

							
Item Size (MM)	DBR47MC 1.2 x 12	DBR48MC 1.8 x 12	DBR49MC 2.5 x 12	DBR50MC 4 x 12	DBR51MC 3 x 1.2	DBR52MC 4 x 1.2	DBR53MC 5 x 1.2




RPM Price 20,000-35,000

			
Item Size (MM)	DBR54MC 3 x 2	DBR55MC 4 x 2	DBR56MC 5 x 2

RPM Price 20,000-35,000

							
Item Size (MM)	DBR57MC 3 x 1.8	DBR58MC 4 x 1.8	DBR59MC 5 x 1.8	DBR60MC 3 x 1	DBR61MC 4 x 1	DBR62MC 5 x 1	DBR63MC 1.2 x 4










RPM Price 20,000-35,000

			
Item Size (MM)	DBR64MC 2.5 x 8	DBR65MC 2.8 x 14	DBR66MC 4.8 x 12

RPM Price 20,000-35,000

Coarse Grit



Fine Grit







									
Item Size (MM)	DBR67C 1.8 x 1.4	DBR68C 2.3 x 1.9	DBR69C 4 x 8	DBR70C 4 x 12	DBR1F 1.8 x 1.4	DBR2F 1.6 x 10	DBR3F 2.5 x 10	DBR4F 1.8 x 8	DBR5F 1.8 x 12





RPM Price 20,000-35,000

# ZIRCONIA / Low Speed Carbide Lab Burs


## Low Speed Carbide Lab Burs

	Round				Inverted Cone			Straight, Flat End, Cross Cut		
Item										
Size (MM)		CBR 05 LAB 006 Ø	CBR 4 LAB 014 Ø	CBR 8 LAB 023 Ø		CBR 34 LAB 008 x 1.0	CBR 35 LAB 010 x 1.0		CBR 556 LAB 009 x 3.8	CBR 557 LAB 010 x 3.8
Price										

	Cross Cut, Super Fine			Cross Cut, Fine		
Item						
Size (MM)	CBR1 SF-CC 16 x 2.3	CBR2 SF-CC 13 x 4.5	CBR3 SF-CC 11.5 x 4	CBR1 F-CC 16 x 2.3	CBR2 F-CC 12 x 6	CBR3 F-CC 14 x 6
Price						

	Cross Cut, Medium								
Item									
Size (MM)	CBR1 M-CC 2.3 x 2.5	CBR2 M-CC 2.8 x 3.0	CBR3 M-CC 3.5 x 4.0	CBR4 M-CC 4.5 x 5.0	CBR5 M-CC 5.5 x 6.0	CBR6 M-CC 7.0 x 8.0	CBR7 M-CC 16.0 x 2.3	CBR8 M-CC 12 x 6	CBR9 M-CC 12 x 6
Price									





























































	Cross Cut, Coarse		Cross Cut, Super Coarse	
Item				
Size (MM)	CBR1 C-CC 12 x 6	CBR2 C-CC 14 x 6	CBR1 SC-CC 12.5 x 6	CBR2 SC-CC 14 x 6
Price				

	Plain Cut Fine		Plain Cut Medium			Plain Cut Coarse	
Item							
Size (MM)	CBR1 F-PC 14 x 4	CBR1 M-PC 2.4 x 2.7	CBR2 M-PC 3 x 4	CBR3 M-PC 5.3 x 6	CBR4 M-PC 12 x 6	CBR5 M-PC 14 x 6	CBR1 C-PC 14 x 6
Price							



**High Speed / Friction Grip Diamond Burs  
Mid-Coarse Grit**





























































**ZIRCONIA / High Speed / Friction Grip Diamond Burs**

										
Item Size (MM)	DBR 1 MCCL 0.8 Ø	DBR 2 MCCL 1 Ø	DBR 3 MCCL 1.2 Ø	DBR 4 MCCL 1.4 Ø	DBR 5 MCCL 1.6 Ø	DBR 6 MCCL 1.8 Ø	DBR 7 MCCL 2 Ø	DBR 8 MCCL 1 Ø	DBR 9 MCCL 1.2 Ø	DBR 10 MCCL 1.4 Ø
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack
										
Item Size (MM)	DBR 11 MCCL 1.2 Ø	DBR 12 MCCL 1.4 Ø	DBR 13 MCCL 1.6 Ø	DBR 14 MCCL 1.2 x 3	DBR 15 MCCL 1.4 x 3.5	DBR 16 MCCL 1.6 x 3.5	DBR 17 MCCL 1.8 x 3.5	DBR 18 MCCL 1.2 x 3	DBR 19 MCCL 1.4 x 3.5	DBR 20 MCCL 1.2 x 10
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack
										
Item Size (MM)	DBR 21 MCCL 1.6 x 10	DBR 22 MCCL 1.8 x 10	DBR 23 MCCL 1.2 x 1.5	DBR 24 MCCL 1.4 x 1.5	DBR 25 MCCL 1.6 x 1.5	DBR 26 MCCL 1.8 x 2.5	DBR 27 MCCL 1.2 x 1.5	DBR 28 MCCL 1.4 x 1.5	DBR 29 MCCL 1.6 x 1.5	DBR 30 MCCL 1.2 x 3.5
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack
										
Item Size (MM)	DBR 31 MCCL 1.4 x 3.5	DBR 32 MCCL 1.6 x 3.5	DBR 33 MCCL 1.8 x 5	DBR 34 MCCL 1.2 x 2.3	DBR 35 MCCL 1.4 x 2.3	DBR 36 MCCL 1.6 x 2.5	DBR 37 MCCL 1.2 x 2.3	DBR 38 MCCL 1.2 x 3.5	DBR 39 MCCL 1 x 3.5	DBR 40 MCCL 1 x 4
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack
										
Item Size (MM)	DBR 41 MCCL 1.4 x 4	DBR 42 MCCL 1.6 x 4	DBR 43 MCCL 1.2 x 6	DBR 44 MCCL 1.4 x 6	DBR 45 MCCL 1.6 x 6	DBR 46 MCCL 1.2 x 8	DBR 47 MCCL 1.4 x 8	DBR 48 MCCL 1.6 x 8	DBR 49 MCCL 1.2 x 7	DBR 50 MCCL 1.4 x 8
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack
										
Item Size (MM)	DBR 51 MCCL 1.6 x 8	DBR 52 MCCL 1.8 x 8	DBR 53 MCCL 1 x 6	DBR 54 MCCL 1.2 x 6	DBR 55 MCCL 1.4 x 6	DBR 56 MCCL 1.2 x 8	DBR 57 MCCL 1.6 x 8	DBR 58 MCCL 1.6 x 10	DBR 59 MCCL 1.8 x 10	DBR 60 MCCL 1.4 x 8
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack

**High Speed / Friction Grip Diamond Burs Continued on Next Page**

**ZIRCONIA / High Speed / Friction Grip Diamond Burs**




















































**Mid-Coarse Grit**

										
Item	DBR 61 MCCL	DBR 62 MCCL	DBR 63 MCCL	DBR 64 MCCL	DBR 65 MCCL	DBR 66 MCCL	DBR 67 MCCL	DBR 68 MCCL	DBR 69 MCCL	DBR 70 MCCL
Size (MM)	2 x 4	2 x 6	2.1 x 10	1 x 4	1.2 x 4	1.4 x 4	1.6 x 4	1.8 x 5	2 x 5	2.3 x 5
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack
										
Item	DBR 71 MCCL	DBR 72 MCCL	DBR 73 MCCL	DBR 74 MCCL	DBR 75 MCCL	DBR 76 MCCL	DBR 77 MCCL	DBR 78 MCCL	DBR 79 MCCL	DBR 80 MCCL
Size (MM)	1.2 x 4	1.4 x 4	1.6 x 4	1.8 x 5	2.3 x 4	1.4 x 6	1.6 x 6	1.4 x 6	1.6 x 6	1 x 8
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack
										
Item	DBR 81 MCCL	DBR 82 MCCL	DBR 83 MCCL	DBR 84 MCCL	DBR 85 MCCL	DBR 86 MCCL	DBR 87 MCCL	DBR 88 MCCL	DBR 89 MCCL	DBR 90 MCCL
Size (MM)	1.2 x 8	1.4 x 8	1.6 x 8	1.8 x 8	2.3 x 9	1 x 8	1.2 x 8	1.4 x 10	1.6 x 10	1.8 x 10
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack
										
Item	DBR 91 MCCL	DBR 92 MCCL	DBR 93 MCCL	DBR 94 MCCL	DBR 95 MCCL	DBR 96 MCCL	DBR 97 MCCL	DBR 98 MCCL	DBR 99 MCCL	DBR 100 MCCL
Size (MM)	1.8 x 12	1 x 5	1.2 x 5	1.4 x 5	1.6 x 5	1.8 x 5	1.6 x 6	1.8 x 6	1.4 x 7	1.6 x 7
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack
										
Item	DBR 101 MCCL	DBR 102 MCCL	DBR 103 MCCL	DBR 104 MCCL	DBR 105 MCCL	DBR 106 MCCL	DBR 107 MCCL	DBR 108 MCCL	DBR 109 MCCL	DBR 110 MCCL
Size (MM)	1.8 x 7	2 x 7	1.6 x 6	1.6 x 7	1.2 x 8	1.4 x 8	1.6 x 8	1.8 x 8	2 x 8	2.4 x 8
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack
										
Item	DBR 111 MCCL	DBR 112 MCCL	DBR 113 MCCL	DBR 114 MCCL	DBR 115 MCCL	DBR 116 MCCL	DBR 117 MCCL	DBR 118 MCCL	DBR 119 MCCL	DBR 120 MCCL
Size (MM)	1.4 x 9	1.6 x 9	1.8 x 9	1.2 x 10	1.4 x 10	1.6 x 10	1.8 x 10	2 x 10	2.2 x 10	1.4 x 11.5
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack











**High Speed / Friction Grip Diamond Burs Continued on Next Page**

**ZIRCONIA / High Speed / Friction Grip Diamond Burs**

**Mid-Coarse Grit**

											
Item Size (MM)	DBR 121 MCCL 2.4 x 11.5	DBR 122 MCCL 1.8 x 10	DBR 123 MCCL 2 x 3.8	DBR 124 MCCL 3 x 5	DBR 125 MCCL 1.2 x 4	DBR 126 MCCL 1.6 x 5	DBR 127 MCCL 1.2 x 7.5	DBR 128 MCCL 1.4 x 7.5	DBR 129 MCCL 1.6 x 7.5	DBR 130 MCCL 1.4 x 7.5	
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	
											
Item Size (MM)	DBR 131 MCCL 1.3 x 10	DBR 132 MCCL 1.6 x 10	DBR 133 MCCL 1.6 x 3.5	DBR 134 MCCL 1.8 x 3.8	DBR 135 MCCL 2.3 x 5	DBR 136 MCCL 2.8 x 5	DBR 137 MCCL 3.2 x 5	DBR 138 MCCL 1.2 x 6	DBR 139 MCCL 1.2 x 6	DBR 140 MCCL 1.2 x 8	
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	
											
Item Size (MM)	DBR 141 MCCL 1.4 x 8	DBR 142 MCCL 1.8 x 6.5	DBR 143 MCCL 1.6 x 8	DBR 144 MCCL 1.8 x 8	DBR 145 MCCL 1.4 x 10	DBR 146 MCCL 1.8 x 10	DBR 147 MCCL 2 x 10	DBR 148 MCCL 1.8 x 1	DBR 149 MCCL 4.5 x 0.6	DBR 150 MCCL 4.0 x 2	
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	
											
Item Size (MM)	DBR 151 MCCL 5 x 2	DBR 152 MCCL 4.2 x 1.7	DBR 153 MCCL 2.3 x 0.8	DBR 154 MCCL 4 x 1.2	DBR 155 MCCL 5 x 1.6	DBR 156 MCCL 1.0 x 4	DBR 157 MCCL 1.2 x 4	DBR 158 MCCL 1.4 x 4	DBR 159 MCCL 1.6 x 5	DBR 160 MCCL 1.8 x 5	
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	
											
Item Size (MM)	DBR 161 MCCL 1.2 x 3	DBR 162 MCCL 2 x 5	DBR 163 MCCL 2.3 x 5	DBR 164 MCCL 3.2 x 5	DBR 165 MCCL 2.3 x 9	DBR 166 MCCL 3.2 x 5	DBR 167 MCCL 3.5 x 7	DBR 168 MCCL 1.2 x 0.5	DBR 169 MCCL 2.9 x 8	DBR 170 MCCL 1.4 x 8	DBR 171 MCCL 1.8 x 10
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack











**Coarse Grit**






										
Item Size (MM)	DBR 1 CCL 1.4 x 1.1	DBR 2 CCL 1.6 x 1.3	DBR 3 CCL 1.8 x 1.4	DBR 4 CCL 1.4 x 1.5	DBR 5 CCL 1.2 x 3.5	DBR 6 CCL 1.0 x 4	DBR 7 CCL 1.4 x 4	DBR 8 CCL 1.2 x 8	DBR 9 CCL 1.4 x 8	DBR 10 CCL 1.6 x 8
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack

High Speed / Friction Grip Diamond Burs Continued on Next Page











**ZIRCONIA / High Speed / Friction Grip Diamond Burs**











**Coarse Grit**











									
Item DBR 11 CCL	Item DBR 12 CCL	Item DBR 13 CCL	Item DBR 14 CCL	Item DBR 15 CCL	Item DBR 16 CCL	Item DBR 17 CCL	Item DBR 18 CCL	Item DBR 19 CCL	Item DBR 20 CCL
Size (MM) 1.6 x 10	Size (MM) 1.6 x 6	Size (MM) 1.6 x 10	Size (MM) 1.8 x 10	Size (MM) 1.8 x 7	Size (MM) 2 x 7	Size (MM) 1.6 x 8	Size (MM) 1.8 x 9	Size (MM) 1.8 x 10	Size (MM) 1.6 x 5
Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack











				
Item DBR 21 CCL	Item DBR 22 CCL	Item DBR 23 CCL	Item DBR 24 CCL	Item DBR 25 CCL
Size (MM) 1.2 x 7.5	Size (MM) 1.8 x 3.8	Size (MM) 2.3 x 5	Size (MM) 4.4 x 1.7	Size (MM) 2 x 9
Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack

**Fine Grit**








									
Item DBR 1 FCL	Item DBR 2 FCL	Item DBR 3 FCL	Item DBR 4 FCL	Item DBR 5 FCL	Item DBR 6 FCL	Item DBR 7 FCL	Item DBR 8 FCL	Item DBR 9 FCL	Item DBR 10 FCL
Size (MM) 0.8 x 0.5	Size (MM) 0.8 x 0.5	Size (MM) 1.2 x 3.5	Size (MM) 1 x 3.5	Size (MM) 0.9 x 2.2	Size (MM) 1.2 x 3.8	Size (MM) 0.9 x 2.2	Size (MM) 1.2 x 8	Size (MM) 1.4 x 8	Size (MM) 1.6 x 10
Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack

									
Item DBR 11 FCL	Item DBR 12 FCL	Item DBR 13 FCL	Item DBR 14 FCL	Item DBR 15 FCL	Item DBR 16 FCL	Item DBR 17 FCL	Item DBR 18 FCL	Item DBR 19 FCL	Item DBR 20 FCL
Size (MM) 1.8 x 10	Size (MM) 1.4 x 6	Size (MM) 1.8 x 8	Size (MM) 1.6 x 10	Size (MM) 1.8 x 10	Size (MM) 1 x 4	Size (MM) 1.4 x 4	Size (MM) 2.1 x 5	Size (MM) 1.6 x 6	Size (MM) 1.8 x 6
Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack
































									
Item DBR 21 FCL	Item DBR 22 FCL	Item DBR 23 FCL	Item DBR 24 FCL	Item DBR 25 FCL	Item DBR 26 FCL	Item DBR 27 FCL	Item DBR 28 FCL	Item DBR 29 FCL	Item DBR 30 FCL
Size (MM) 1.4 x 6	Size (MM) 1.4 x 7	Size (MM) 1.6 x 7	Size (MM) 1.8 x 7	Size (MM) 1.2 x 8	Size (MM) 1.4 x 8	Size (MM) 1.6 x 8	Size (MM) 1.8 x 8	Size (MM) 1.8 x 9	Size (MM) 1.4 x 10
Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack

									
Item DBR 31 FCL	Item DBR 32 FCL	Item DBR 33 FCL	Item DBR 34 FCL	Item DBR 35 FCL	Item DBR 36 FCL	Item DBR 37 FCL	Item DBR 38 FCL	Item DBR 39 FCL	Item DBR 40 FCL
Size (MM) 1.6 x 10	Size (MM) 1.8 x 10	Size (MM) 2 x 10	Size (MM) 0.9 x 4	Size (MM) 1.2 x 7.5	Size (MM) 1.4 x 7.5	Size (MM) 1.6 x 7.5	Size (MM) 1.6 x 3.5	Size (MM) 1.8 x 3.8	Size (MM) 2.3 x 5
Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack	Price 5-Pack

**Fine Grit**









































							
Item	DBR 41 FCL	DBR 42 FCL	DBR 43 FCL	DBR 44 FCL	DBR 45 FCL	DBR 46 FCL	DBR 47 FCL
Size (MM)	1.4 x 8	1.6 x 8	1.8 x 8	1.8 x 0.8	2 x 5	2.3 x 5	1.2 x 0.5
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack

**Extra-Fine Grit**

											
Item	DBR 48 EFCL	DBR 49 EFCL	DBR 50 EFCL	DBR 51 EFCL	DBR 52 EFCL	DBR 53 EFCL	DBR 54 EFCL	DBR 55 EFCL	DBR 56 EFCL	DBR 57 EFCL	
Size (MM)	1 x 0.7	1 x 6	1.2 x 6	1.6 x 6	1.6 x 10	1.4 x 5	2.1 x 5	1.6 x 3.5	1.8 x 6	1.6 x 7	
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	
											
Item	DBR 58 EFCL	DBR 59 EFCL	DBR 60 EFCL	DBR 61 EFCL	DBR 62 EFCL	DBR 63 EFCL	DBR 64 EFCL	DBR 65 EFCL	DBR 66 EFCL	DBR 67 EFCL	
Size (MM)	1.2 x 8	1.6 x 8	1.6 x 9	1.6 x 10	1.8 x 10	2 x 10	0.9 x 3.5	1.2 x 4	1.2 x 7.5	1.4 x 7.5	
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	
											
Item	DBR 68 EFCL	DBR 69 EFCL	DBR 70 EFCL	DBR 71 EFCL	DBR 72 EFCL	DBR 73 EFCL	DBR 74 EFCL	DBR 75 EFCL	DBR 76 EFCL	DBR 77 EFCL	DBR 78 EFCL
Size (MM)	1.6 x 7.5	1.3 x 10	1.6 x 10	1.6 x 3.5	1.8 x 3.8	2.3 x 5	1.6 x 8	1.8 x 8	2 x 5	2.3 x 5	3.2 x 5
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack

# ZIRCONIA / High Speed / Friction Grip Carbide Burs

## High Speed / Friction Grip Carbide Burs

















										
Item	CBR05	CBR1	CBR2	CBR4	CBR6	CBR8		CBR34	CBR35	CBR38
Size (MM)	005 Ø	008 Ø	010 Ø	014 Ø	018 Ø	023 Ø		008 x 0.8	010 x 1.0	018 x 1.8
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack		5-Pack	5-Pack	5-Pack
										
Item	CBR56	CBR57		CBR169	CBR169L	CBR170	CBR170L		CBR330	
Size (MM)	009 x 4.0	010 x 4.5		009 x 4.0	009 x 5.0	010 x 4.5	010 x 5.5?		008 x 1.7	
Price	5-Pack	5-Pack		5-Pack	5-Pack	5-Pack	5-Pack		5-Pack	
										
Item	CBR556	CBR557	CBR557L	CBR558		CBR699	CBR699L	CBR700	CBR700L	CBR701
Size (MM)	009 x 4.0	010 x 4.5	010 x 6.0	012 x 4.5		009 x 4.0	009 x 5.0	010 x 4.5	010 x 6.0	012 x 4.5
Price	5-Pack	5-Pack	5-Pack	5-Pack		5-Pack	5-Pack	5-Pack	5-Pack	5-Pack
										
Item	CBR1156		CBR1171		CBR7204F	CBR7205F	CBR7206F		CBR7714F	
Size (MM)	009 x 4.0		012 x 5.0		014 x 9.5	016 x 9.5	018 x 9.5		014 x 6.0	
Price	5-Pack		5-Pack		5-Pack	5-Pack	5-Pack		5-Pack	
										
Item	CBR7801F	CBR0151DC	CBR0152DC							
Size (MM)	009 x 4.0	015 x 11.0	015 x 9.0							
Price	5-Pack	5-Pack	5-Pack							

**SPC: Silicone Polishers for Ceramic and Resin  
Coarse Grit**

**CERAMICS / Silicone Polishers**

								
Item	SPC1C	SPC2C	SPC3C	SPC4C	SPC5C	SPC6C	SPC7C	SPC8C
Size (MM)	16 x 1	16 x 2	16 x 2	22 x 1	22 x 3	22 x 3.5	6 x 23	7 x 20
MAX RPM	15,000	15,000	15,000	10,000	10,000	10,000	20,000	20,000
Price								
								
Item	SPC9C	SPC10C	SPC11C	SPC12C	SPC13C	SPC14C	SPC15C	SPC16C
Size (MM)	3 x 7	5 x 2	3 x 7	4 x 12	6 x 16	4 x 12	9 x 8.5	18 x 3
MAX RPM	20,000	20,000	20,000	20,000	20,000	20,000	12,000	12,000
Price								
								
Item	SPC17C	SPC18C	SPC19C	SPC20C	SPC21C	SPC22C	SPC23C	
Size (MM)	12 x 1	18 x 1	12 x 2	15 x 2	12 x 4	12 x 2	15 x 2	
MAX RPM	20,000	20,000	20,000	20,000	20,000	20,000	20,000	
Price								

**Mid-Coarse Grit**

								
Item	SPC1MC	SPC2MC	SPC3MC	SPC4MC	SPC5MC	SPC6MC	SPC7MC	SPC8MC
Size (MM)	16 x 1	16 x 2	16 x 2	22 x 1	22 x 3	22 x 3.5	6 x 23	7 x 20
MAX RPM	15,000	15,000	15,000	10,000	10,000	10,000	20,000	20,000
Price								
								
Item	SPC9MC	SPC10MC	SPC11MC	SPC12MC	SPC13MC	SPC14MC	SPC15MC	SPC16MC
Size (MM)	3 x 7	5 x 2	3 x 7	4 x 12	6 x 16	4 x 12	9 x 8.5	18 x 3
MAX RPM	15,000	15,000	15,000	15,000	15,000	15,000	10,000	10,000
Price								









Silicone Polishers for Ceramic and Resin Continued on Next Page









# CERAMICS / Silicone Polishers

## Mid-Coarse Grit

							
Item	SPC17MC	SPC18MC	SPC19MC	SPC20MC	SPC21MC	SPC22MC	SPC23MC
Size (MM)	12 x 1	18 x 1	12 x 2	15 x 2	12 x 4	12 x 2	15 x 2
MAX RPM	15,000	15,000	15,000	15,000	15,000	15,000	15,000
Price							

## Fine Grit

								
Item	SPC1F	SPC2F	SPC3F	SPC4F	SPC5F	SPC6F	SPC7F	SPC8F
Size (MM)	16 x 1	16 x 2	16 x 2	22 x 1	22 x 3	22 x 3.5	6 x 23	7 x 20
MAX RPM	10,000	10,000	10,000	8,000	8,000	8,000	10,000	10,000
Price								









								
Item	SPC9F	SPC10F	SPC11F	SPC12F	SPC13F	SPC14F	SPC15F	SPC16F
Size (MM)	3 x 7	5 x 2	3 x 7	4 x 12	6 x 16	4 x 12	9 x 8.5	18 x 3
MAX RPM	8,000	8,000	8,000	8,000	8,000	8,000	8,000	8,000
Price								

							
Item	SPC17F	SPC18F	SPC19F	SPC20F	SPC21F	SPC22F	SPC23F
Size (MM)	12 x 1	18 x 1	12 x 2	15 x 2	12 x 4	12 x 2	15 x 2
MAX RPM	8,000	8,000	8,000	8,000	8,000	8,000	8,000
Price							



# ALLOYS / Silicone Polishers









## SPA: Silicone Polishers for Alloy Coarse Grit








								
Item	SPA1C	SPA2C	SPA3C	SPA4C	SPA5C	SPA6C	SPA7C	SPA8C
Size (MM)	16 x 1	16 x 2	16 x 2	22 x 1	22 x 3	22 x 3.5	6 x 23	7 x 20
MAX RPM	8,000-10,000							
Price								

								
Item	SPA9C	SPA10C	SPA11C	SPA12C	SPA13C	SPA14C	SPA15C	SPA16C
Size (MM)	3 x 7	5 x 12	3 x 7	4 x 12	6 x 16	4 x 12	9 x 8.5	18 x 3
MAX RPM	10,000-12,000							
Price								

							
Item	SPA17C	SPA18C	SPA19C	SPA20C	SPA21C	SPA22C	SPA23C
Size (MM)	12 x 1	18 x 1	12 x 2	15 x 2	12 x 4	12 x 2	15 x 2
MAX RPM	10,000-12,000						
Price							

## Mid-Coarse Grit

								
Item	SPA1MC	SPA2MC	SPA3MC	SPA4MC	SPA5MC	SPA6MC	SPA7MC	SPA8MC
Size (MM)	16 x 1	16 x 2	16 x 2	22 x 1	22 x 3	22 x 3.5	6 x 23	7 x 20
MAX RPM	8,000-10,000							
Price								

								
Item	SPA9MC	SPA10MC	SPA11MC	SPA12MC	SPA13MC	SPA14MC	SPA15MC	SPA16MC
Size (MM)	3 x 7	5 x 12	3 x 7	4 x 12	6 x 16	4 x 12	9 x 8.5	18 x 3
MAX RPM	10,000-12,000							
Price								

# ALLOYS / Silicone Polishers

## Mid-Coarse Grit



Item	SPA17MC	SPA18MC	SPA19MC	SPA20MC	SPA21MC	SPA22MC	SPA23MC
Size (MM)	12 x 1	18 x 1	12 x 2	15 x 2	12 x 4	12 x 2	15 x 2
MAX RPM	10,000-12,000						
Price							

## Fine Grit



Item	SPA1F	SPA2F	SPA3F	SPA4F	SPA5F	SPA6F	SPA7F	SPA8F
Size (MM)	16 x 1	16 x 2	16 x 2	22 x 1	22 x 3	22 x 3.5	6 x 23	7 x 20
MAX RPM	8,000-10,000							
Price								



Item	SPA9F	SPA10F	SPA11F	SPA12F	SPA13F	SPA14F	SPA15F	SPA16F
Size (MM)	3 x 7	5 x 12	3 x 7	4 x 12	6 x 16	4 x 12	9 x 8.5	18 x 3
MAX RPM	10,000-12,000							
Price								



Item	SPA17F	SPA18F	SPA19F	SPA20F	SPA21F	SPA22F	SPA23F
Size (MM)	12 x 1	18 x 1	12 x 2	15 x 2	12 x 4	12 x 2	15 x 2
MAX RPM	10,000-12,000						
Price							
















## Extra Fine Grit


























Item	SPA1EF	SPA2EF	SPA3EF	SPA4EF	SPA5EF	SPA6EF	SPA7EF	SPA8EF
Size (MM)	16 x 1	16 x 2	16 x 2	22 x 1	22 x 3	22 x 3.5	6 x 23	7 x 20
MAX RPM	8,000-10,000							
Price								

**ALLOYS / Silicone  
Polishers**

**Extra Fine Grit**

								
Item	SPA9EF	SPA10EF	SPA11EF	SPA12EF	SPA13EF	SPA14EF	SPA15EF	SPA16EF
Size (MM)	3 x 7	5 x 12	3 x 7	4 x 12	6 x 16	4 x 12	9 x 8.5	18 x 3
MAX RPM	10,000-12,000							
Price								
								
Item	SPA17EF	SPA18EF	SPA19EF	SPA20EF	SPA21EF	SPA22EF	SPA23EF	
Size (MM)	12 x 1	18 x 1	12 x 2	15 x 2	12 x 4	12 x 2	15 x 2	
MAX RPM	10,000-12,000							
Price								

**RPA: Rubber Polishers for Alloy / Extra Coarse Grit**

								
Item	RPA1C	RPA2C	RPA3C	RPA4C	RPA5C	RPA6C	RPA7C	RPA8C
Size (MM)	16 x 1	16 x 2	16 x 2	22 x 1	22 x 3	22 x 3.5	6 x 23	7 x 20
MAX RPM	15,000	15,000	15,000	10,000	10,000	10,000	20,000	20,000
Price								
								
Item	RPA1MC	RPA2MC	RPA3MC	RPA4MC	RPA5MC	RPA6MC	RPA7MC	RPA8MC
Size (MM)	16 x 1	16 x 2	16 x 2	22 x 1	22 x 3	22 x 3.5	6 x 23	7 x 20
MAX RPM	12,000	12,000	12,000	8,000	8,000	8,000	15,000	15,000
Price								
								
Item	RPA2EC	RPA3EC	RPA4EC	RPA5EC	RPA6EC	RPA7EC	RPA8EC	
Size (MM)	16 x 2	16 x 2	22 x 1	22 x 3	22 x 3.5	6 x 23	7 x 20	
MAX RPM	15,000	15,000	10,000	10,000	10,000	20,000	20,000	
Price								

**Rubber Polishers for Alloy Continued on Next Page**

*OpaDental.com Customer Service: 304-523-6721*

# ALLOYS / Rubber Polishers

## Fine Grit

Item	RPA1F	RPA2F	RPA3F	RPA4F	RPA5F	RPA6F	RPA7F	RPA8F
Size (MM)	16 x 1	16 x 2	16 x 2	22 x 1	22 x 3	22 x 3.5	6 x 23	7 x 20
MAX RPM	10,000	10,000	10,000	8,000	8,000	8,000	10,000	10,000
Price								

## Coarse Grit

Item	RPA9C	RPA10C	RPA11C	RPA12C	RPA13C	RPA14C	RPA15C	RPA16C	RPA17C
Size (MM)	3 x 7	5 x 12	5 x 16	3 x 7	4 x 12	6 x 16	6 x 16	9 x 8.5	18 x 3
MAX RPM	20,000	20,000	20,000	20,000	20,000	20,000	20,000	12,000	12,000
Price									

Item	RPA18C	RPA19C	RPA20C	RPA21C	RPA22C	RPA23C	RPA24C	RPA25C
Size (MM)	12 x 1	18 x 1	23 x 1	12 x 2	15 x 2	12 x 4	12 x 2	15 x 2
MAX RPM	20,000	20,000	20,000	20,000	20,000	20,000	20,000	20,000
Price								

## Fine Grit










Item	RPA9F	RPA10F	RPA11F	RPA12F	RPA13F	RPA14F	RPA15F	RPA16F	RPA17F
Size (MM)	3 x 7	5 x 12	5 x 16	3 x 7	4 x 12	6 x 16	6 x 16	9 x 8.5	18 x 3
MAX RPM	15,000	15,000	15,000	15,000	15,000	15,000	15,000	12,000	12,000
Price									

Item	RPA18F	RPA19F	RPA20F	RPA21F	RPA22F	RPA23F	RPA24F	RPA25F
Size (MM)	12 x 1	18 x 1	23 x 1	12 x 2	15 x 2	12 x 4	12 x 2	15 x 2
MAX RPM	15,000	15,000	15,000	15,000	15,000	15,000	15,000	15,000
Price								

Rubber Polishers for Alloy continued on Next Page










# ALLOYS / Rubber Polishers



## Extra Fine Grit

									
Item	RPA1EF	RPA2EF	RPA3EF	RPA4EF	RPA5EF	RPA6EF	RPA7EF	RPA8EF	RPA9EF
Size (MM)	3 x 7	5 x 12	5 x 16	3 x 7	4 x 12	6 x 16	6 x 16	9 x 8.5	18 x 3
MAX RPM	12,000	12,000	12,000	12,000	12,000	12,000	12,000	10,000	10,000
Price									

								
Item	RPA10EF	RPA11EF	RPA12EF	RPA13EF	RPA14EF	RPA15EF	RPA16EF	RPA17EF
Size (MM)	12 x 1	18 x 1	23 x 1	12 x 2	15 x 2	12 x 4	12 x 2	15 x 2
MAX RPM	12,000	12,000	12,000	12,000	12,000	12,000	12,000	12,000
Price								

## Ultra Fine Grit

									
Item	RPA1UF	RPA2UF	RPA3UF	RPA4UF	RPA5UF	RPA6UF	RPA7UF	RPA8UF	RPA9UF
Size (MM)	3 x 7	5 x 12	5 x 16	3 x 7	4 x 12	6 x 16	6 x 16	9 x 8.5	18 x 3
MAX RPM	8,000	8,000	8,000	8,000	8,000	8,000	8,000	8,000	8,000
Price									

								
Item	RPA10UF	RPA11UF	RPA12UF	RPA13UF	RPA14UF	RPA15UF	RPA16UF	RPA17UF
Size (MM)	12 x 1	18 x 1	23 x 1	12 x 2	15 x 2	12 x 4	12 x 2	15 x 2
MAX RPM	8,000	8,000	8,000	8,000	8,000	8,000	8,000	8,000
Price								

# ACRYLIC / Silicone Polishers

## SPR: Silicone Polishers for Resin / Acrylic Extra Coarse Grit

								
Item Size (MM)	SPR1EC	SPR2EC	SPR3EC	SPR4EC	SPR5EC	SPR6EC	SPR7EC	SPR8EC
MAX RPM	4 x 12	6 x 16	10 x 20	10 x 24	10 x 17	13 x 20	14 x 17	15 x 17
Price	20,000	20,000	12,000	12,000	12,000	10,000	10,000	10,000

## Coarse Grit

								
Item Size (MM)	SPR9C	SPR10C	SPR11C	SPR12C	SPR13C	SPR14C	SPR15C	SPR16C
MAX RPM	4 x 12	6 x 16	10 x 20	10 x 24	10 x 17	13 x 20	14 x 17	15 x 17
Price	20,000	20,000	12,000	12,000	12,000	10,000	10,000	10,000

## Mid-Coarse Grit


								
Item Size (MM)	SPR17MC	SPR18MC	SPR19MC	SPR20MC	SPR21MC	SPR22MC	SPR23MC	SPR24MC
MAX RPM	4 x 12	6 x 16	10 x 20	10 x 24	10 x 17	13 x 20	14 x 17	15 x 17
Price	15,000	15,000	10,000	10,000	10,000	8,000	8,000	8,000

## Extra Fine Grit

								
Item Size (MM)	SPR25EF	SPR26EF	SPR27EF	SPR28EF	SPR29EF	SPR30EF	SPR31EF	SPR32EF
MAX RPM	4 x 12	6 x 16	10 x 20	10 x 24	10 x 17	13 x 20	14 x 17	15 x 17
Price	7,000	7,000	7,000	7,000	7,000	6,000	6,000	6,000

# CLINICAL PRODUCTS

## PKDZL: Clinical Chairside Kits



Coarse

Fine

The Chairside Kits are organized so that you start from **left to right**, depending on the extent of the adjustments you need to perform. We recommend that you begin with the coarse instruments first for two reasons: your rotary instruments will last longer, and you will have to work less because you are letting the instrument do the work for you!


Chairside Kits are available for the following materials:

- Zirconia (Traditional and Cosmetic)
- Lithium Disilicate (e.max)
- Porcelain
- Metal
- Gold
- Composite

### Customizable:



You can make your own custom kit by substituting different instruments of the same type, e.g. diamond bur “22MCCL” for diamond bur “34MCCL”, etc. You can also replace individual pieces as they wear out instead of buying a new kit!

**Clinical Chairside Kits Continued on Next Page**



**Item: PKDZL  
Chairside Kit For Zirconia Contents:**

**Ceramic Diamond Grinders for Zirconia, Latch**



		
Item:	<b>CDZ 2 MCCL</b>	<b>CDZ 4 MCCL</b>
Size:	<b>3.6 x 7</b>	<b>4 x 10.5</b>
Qty:	<b>1</b>	<b>1</b>

**Rubber Diamond Polishers for Zirconia, Latch**








**Mid-Coarse Grit**

			
Item:	<b>RDZ 14 MCCL</b>	<b>RDZ 18 MCCL</b>	<b>RDZ 19 MCCL</b>
Size:	<b>3.5 x 10.5</b>	<b>10 x 2</b>	<b>10 x 1.8</b>
Qty:	<b>1</b>	<b>1</b>	<b>1</b>

**Extra Fine Grit**

			
Item:	<b>RDZ 32 EFCL</b>	<b>RDZ 38 EFCL</b>	<b>RDZ 39 EFCL</b>
Size:	<b>4 x 12</b>	<b>10 x 2</b>	<b>10 x 1.8</b>
Qty:	<b>1</b>	<b>1</b>	<b>1</b>






**High Speed / Friction Grip Diamond Burs**

							
Item:	<b>DBR 3 CCL</b>	<b>DBR 24 MCCL</b>	<b>DBR 113 MCCL</b>	<b>DBR 124 MCCL</b>	<b>DBR 148 MCCL</b>	<b>DBR 154 MCCL</b>	<b>DBR 40 FCL</b>
Size:	<b>1.8 Ø</b>	<b>1.4 x 1.5</b>	<b>1.8 x 9</b>	<b>3 x 5</b>	<b>1.8 x 1</b>	<b>4 x 1.2</b>	<b>2.3 x 5</b>
Qty:	<b>1</b>	<b>1</b>	<b>1</b>	<b>1</b>	<b>1</b>	<b>1</b>	<b>1</b>







**CLINICAL / Ceramic Diamond Grinders**

**Ceramic Diamond Grinders for Zirconia (CDZ)  
Mid-Coarse Grit**

						
Item	CDZ 1 MCCL	CDZ 2 MCCL	CDZ 3 MCCL	CDZ 4 MCCL	CDZ 5 MCCL	CDZ 6 MCCL
Size (MM)	Ø 3	3.6 x 7	3 x 7	4 x 10.5	3 x 7	3 x 7
RPM	8,000-12,000					
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack


**Ceramic Grinders for Ceramic / Resin / Alloy  
Mid-Coarse Grit**

				
Item	CCA 1 MCCL	CCA 2 MCCL	CCA 3 MCCL	CCA 4 MCCL
Size (MM)	3 x 3	2.5 x 6	2.8 x 7	3 x 7
RPM	8,000-12,000			
Price	5-Pack	5-Pack	5-Pack	5-Pack

**Mid-Coarse Grit, Latch**

								
Item	CCA 5 MCCL	CCA 6 MCCL	CCA 7 MCCL	CCA 8 MCCL	CCA 9 MCCL	CCA 10 MCCL	CCA 11 MCCL	CCA 12 MCCL
Size (MM)	3 x 3	8 x 2	6 x 1.5	3.6 x 7	2.5 x 6	2.8 x 7	3 x 7	2.5 x 6
RPM	8,000-12,000							
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack

**Fine Grit (Non-Precious Alloy)**

				
Item	ACA 1 FCL	ACA 2 FCL	ACA 3 FCL	ACA 4 FCL
Size (MM)	3 x 3	2.5 x 6	2.8 x 7	3 x 7
RPM	8,000-12,000			
Price	5-Pack	5-Pack	5-Pack	5-Pack

**Fine Grit, Latch (Non-Precious Alloy)**

								
Item	ACA 5 FCL	ACA 6 FCL	ACA 7 FCL	ACA 8 FCL	ACA 9 FCL	ACA 10 FCL	ACA 11 FCL	ACA 12 FCL
Size (MM)	3 x 3	8 x 2	6 x 1.5	3.6 x 7	2.5 x 6	2.8 x 7	3 x 7	2.5 x 6
RPM	8,000-12,000							
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack

**CLINICAL / Ceramic Grinders for Alloy**

# CLINICAL / Rubber Diamond Polishers for Zirconia

## Clinical Rubber Diamond Polishers for Zirconia

### Mid-Coarse Grit

### Fine Grit

### Extra Fine Grit

									
Item	RDZ 1 MCCL	RDZ 2 MCCL	RDZ 3 MCCL	RDZ 4 FCL	RDZ 5 FCL	RDZ 6 FCL	RDZ 7 EFCL	RDZ 8 EFCL	RDZ 9 EFCL
Size (MM)	3 x 7	3 x 7	5 x 10	3 x 7	3 x 7	5 x 10	3 x 7	3 x 7	5 x 10
RPM	7,000-10,000								
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack

### Mid-Coarse Grit, Latch

										
Item	RDZ 10 MCCL	RDZ 11 MCCL	RDZ 12 MCCL	RDZ 13 MCCL	RDZ 14 MCCL	RDZ 15 MCCL	RDZ 16 MCCL	RDZ 17 MCCL	RDZ 18 MCCL	RDZ 19 MCCL
Size (MM)	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9 x 8.5	10 x 1	10 x 2	10 x 1.8
RPM	7,000-10,000									
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack

### Fine Grit, Latch

										
Item	RDZ 20 FCL	RDZ 21 FCL	RDZ 22 FCL	RDZ 23 FCL	RDZ 24 FCL	RDZ 25 FCL	RDZ 26 FCL	RDZ 27 FCL	RDZ 28 FCL	RDZ 29 FCL
Size (MM)	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9 x 8.5	10 x 1	10 x 2	10 x 1.8
RPM	7,000-10,000									
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack

### Extra Fine Grit, Latch

										
Item	RDZ 30 EFCL	RDZ 31 EFCL	RDZ 32 EFCL	RDZ 33 EFCL	RDZ 34 EFCL	RDZ 35 EFCL	RDZ 36 EFCL	RDZ 37 EFCL	RDZ 38 EFCL	RDZ 39 EFCL
Size (MM)	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9 x 8.5	10 x 1	10 x 2	10 x 1.8
RPM	7,000-10,000									
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack




# CLINICAL / Rubber Polishers for Alloy

## Clinical Rubber Polishers for Alloy

	Coarse Grit			Fine Grit		Extra Fine Grit			Ultra Fine Grit			
Item	RA 1 CCL	RA 2 CCL	RA 3 CCL	RA 4 FCL	RA 5 FCL	RA 6 FCL	RA 7 EFCL	RA 8 EFCL	RA 9 EFCL	RA 10 UFCL	RA 11 UFCL	RA 12 UFCL
Size (MM)	3 x 7	3 x 7	5 x 10	3 x 7	3 x 7	5 x 10	3 x 7	3 x 7	5 x 10	3 x 7	3 x 7	5 x 10
RPM						7,000-10,000						
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack
<b>Coarse Grit, Latch</b>												
Item	RA 13 CCL	RA 14 CCL	RA 15 CCL	RA 16 CCL	RA 17 CCL	RA 18 CCL	RA 19 CCL	RA 20 CCL	RA 21 CCL	RA 22 CCL		
Size (MM)	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9 x 8.5	10 x 1	10 x 2	10 x 1.8		
RPM					7,000-10,000							
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack		
<b>Fine Grit, Latch</b>												
Item	RA 23 FCL	RA 24 FCL	RA 25 FCL	RA 26 FCL	RA 27 FCL	RA 28 FCL	RA 29 FCL	RA 30 FCL	RA 31 FCL	RA 32 FCL		
Size (MM)	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9 x 8.5	10 x 1	10 x 2	10 x 1.8		
RPM					7,000-10,000							
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack		
<b>Extra Fine Grit, Latch</b>												
Item	RA 33 EFCL	RA 34 EFCL	RA 35 EFCL	RA 36 EFCL	RA 37 EFCL	RA 38 EFCL	RA 39 EFCL	RA 40 EFCL	RA 41 EFCL	RA 42 EFCL		
Size (MM)	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9 x 8.5	10 x 1	10 x 2	10 x 1.8		
RPM					7,000-10,000							
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack		

# CLINICAL / Rubber Polishers for Alloy

## Ultra Fine Grit, Latch

										
Item	RA 43 UFCL	RA 44 UFCL	RA 45 UFCL	RA 46 UFCL	RA 47 UFCL	RA 48 UFCL	RA 49 UFCL	RA 50 UFCL	RA 51 UFCL	RA 52 UFCL
Size (MM)	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9 x 8.5	10 x 1	10 x 2	10 x 1.8
RPM					7,000-10,000					
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack

**CLINICAL / Silicone Diamond Polishers for Resin**

**Clinical Silicone Diamond Polishers for Resin**

**Extra Fine Grit**

**Ultra Fine Grit**

						
Item	SDR 1 EFCL	SDR 2 EFCL	SDR 3 EFCL	SDR 14 UFCL	SDR 15 UFCL	SDR 16 UFCL
Size (MM)	3 x 7	3 x 7	5 x 10	3 x 7	3 x 7	5 x 10
RPM				7,000-10,000		
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack

**Extra Fine Grit, Latch**

										
Item	SDR 4 EFCL	SDR 5 EFCL	SDR 6 EFCL	SDR 7 EFCL	SDR 8 EFCL	SDR 9 EFCL	SDR 10 EFCL	SDR 11 EFCL	SDR 12 EFCL	SDR 13 EFCL
Size (MM)	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9 x 8.5	10 x 1	10 x 2	10 x 1.8
RPM					7,000-10,000					
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack

**Ultra Fine Grit, Latch**

										
Item	SDR 17 UFCL	SDR 18 UFCL	SDR 19 UFCL	SDR 20 UFCL	SDR 21 UFCL	SDR 22 UFCL	SDR 23 UFCL	SDR 24 UFCL	SDR 25 UFCL	SDR 26 UFCL
Size (MM)	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9 x 8.5	10 x 1	10 x 2	10 x 1.8
RPM					7,000-10,000					
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack

**Clinical Silicone Polishers for Ceramic & Resin**

**Coarse Grit**

**Mid-Coarse Grit**

**Fine Grit**

									
Item	SCR 1 CCL	SCR 2 CCL	SCR 3 CCL	SCR 4 MCCL	SCR 5 MCCL	SCR 6 MCCL	SCR 7 FCL	SCR 8 FCL	SCR 9 FCL
Size (MM)	3 x 7	3 x 7	5 x 10	3 x 7	3 x 7	5 x 10	3 x 7	3 x 7	5 x 10
RPM					7,000-10,000				
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack

**Coarse Grit, Latch**

										
Item	SCR 10 CCL	SCR 11 CCL	SCR 12 CCL	SCR 13 CCL	SCR 14 CCL	SCR 15 CCL	SCR 16 CCL	SCR 17 CCL	SCR 18 CCL	SCR 19 CCL
Size (MM)	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9 x 8.5	10 x 1	10 x 2	10 x 1.8
RPM					7,000-10,000					
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack

**Mid-Coarse Grit, Latch**

										
Item	SCR 20 MCCL	SCR 21 MCCL	SCR 22 MCCL	SCR 23 MCCL	SCR 24 MCCL	SCR 25 MCCL	SCR 26 MCCL	SCR 27 MCCL	SCR 28 MCCL	SCR 29 MCCL
Size (MM)	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9 x 8.5	10 x 1	10 x 2	10 x 1.8
RPM					7,000-10,000					
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack

**Fine Grit, Latch**

										
Item	SCR 30 FCL	SCR 31 FCL	SCR 32 FCL	SCR 33 FCL	SCR 34 FCL	SCR 35 FCL	SCR 36 FCL	SCR 37 FCL	SCR 38 FCL	SCR 39 FCL
Size (MM)	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9 x 8.5	10 x 1	10 x 2	10 x 1.8
RPM					7,000-10,000					
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack

**Clinical / Silicone Polishers for Noble / High Noble Alloy**

**Clinical Silicone Polishers for Noble/High Noble Alloy**

**Coarse Grit**

									
Item	SA 1 CCL	SA 2 CCL	SA 3 CCL	SA 4 MCCL	SA 5 MCCL	SA 6 MCCL	SA 7 FCL	SA 8 FCL	SA 9 FCL
Size (MM)	3 x 7	3 x 7	5 x 10	3 x 7	3 x 7	5 x 10	3 x 7	3 x 7	5 x 10
RPM					7,000-10,000				
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack

**Coarse Grit, Latch**

										
Item	SA 10 CCL	SA 11 CCL	SA 12 CCL	SA 13 CCL	SA 14 CCL	SA 15 CCL	SA 16 CCL	SA 17 CCL	SA 18 CCL	SA 19 CCL
Size (MM)	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9 x 8.5	10 x 1	10 x 2	10 x 1.8
RPM					7,000-10,000					
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack

**Mid-Coarse Grit, Latch**

										
Item	SA 20 MCCL	SA 21 MCCL	SA 22 MCCL	SA 23 MCCL	SA 24 MCCL	SA 25 MCCL	SA 26 MCCL	SA 27 MCCL	SA 28 MCCL	SA 29 MCCL
Size (MM)	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9 x 8.5	10 x 1	10 x 2	10 x 1.8
RPM					7,000-10,000					
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack


**Fine Grit, Latch**

										
Item	SA 30 FCL	SA 31 FCL	SA 32 FCL	SA 33 FCL	SA 34 FCL	SA 35 FCL	SA 36 FCL	SA 37 FCL	SA 38 FCL	SA 39 FCL
Size (MM)	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9 x 8.5	10 x 1	10 x 2	10 x 1.8
RPM					7,000-10,000					
Price	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack	5-Pack

# **OpaScan** **Digital Dental System**

## What can you scan?

- ⌘ Impressions
- ⌘ Models
- ⌘ Implants
- ⌘ Single Dies
- ⌘ Multi-Dies
- ⌘ Crown and Bridge
- ⌘ Articulators
- ⌘ Full Arch
- ⌘ Texture Scan



## What will your scan create?

- ⌘ Crowns
- ⌘ Bridges
- ⌘ Veneers
- ⌘ Implant Abutments
- ⌘ Tissue Shaded Custom Healing Abutments
- ⌘ Surgical Stents
- ⌘ Temporaries
- ⌘ Partials
- ⌘ Dentures

## Features

Compact, fast, and efficient, the OpaScan compliments any operation. Utilizing high precision blue and red light technology, the OpaScan features automatic calibration, alignment, and single click scanning.

The OpaScan is remarkably user-friendly, and will make a great addition to any dental organization.



# OpaFit

Dental Impression Materials



## OpaFit 1P

**Putty Base (400g) + Catalyst (400g)**  
Includes: One 400g Base jar, One 400g Catalyst jar, 2 mixing spoons



## OpaFit 2LC

**Light Body: Double Barrel Bio-Cartridges**  
Includes: Two 50ml Cartridges, 6 mixers, 6 tips



## OpaFit 3HC

**Heavy Body: Double Barrel Bio-Cartridges**  
Includes: Two 50ml Cartridges, 6 mixers

**OpaFit Accessories:** Dispenser Gun, Additional Mixers, Additional Tips available